




HiToP[®] Therapy

Indications with HiToP[®]



Index

Origin of the HighTone Power Therapy	5
HiToP® programs	7
Dermatology	10
Neurodermatitis	10
Psoriasis	11
Radiation injuries	12
Gynecology and urology	13
Adnexitis, chronic	13
Cystitis, Cystopyelitis, Cystopyelonephritis, chronische	14
Dysmenorrhea	15
Pain caused by IUPs	16
Stress-urine incontinence	17
Otorhinolaryngology	18
Laryngitis , chronic	18
Chronic otitis media	19
Frontal sinusitis , chronic	20
Maxillary sinusitis , chronic	21
Tinnitus	22
Internal medicin and angiology	23
Adipositas, overweight	23
Arterial hypertonia Arterial occlusive diseases	24
Uratic arthritis	25
Spasmodic asthma	26
Decubitus in the back region	27
Decubitus in the region of the heel bone	28
Diabetic angiopathy	29
Thromboangiitis obliterans	30
Hematoma in the sural region	31
Cardiac discomforts, functional	32
Migraine (treatment at intervals)	33
Migraine and BING-HORTON syndrome	34
Obstipation atonic Obstipation, spastic	35
Posttraumatic edema in the area of the hand	36
Posttraumatic edema in the area of the foot	37
Raynaud's syndrome	38
Thrombosis prophylaxis	39
Stasis ulcer	40
Venous insufficiency , chronic	41
Neurology	42
Amyotrophic lateral sklerosis	42
Atypical faciocephalalgia	43
Diabetic polyneuropathy	44
Nocturnal enuresis	45
Facial nerve paresis	46
Polyneuropathies	47
Shoulder pain in hemiplegia	48
Tension headache, Vasomotoric headache	49
Ophthalmology	50
Conjunctivitis, allergic	50
Orthopaedics, surgery, sports medicine	51
Arthrosis of the carpus	51
Arthrosis of the hip joints, Coxarthrosis	52
Arthrosis of the knee joints	53
Arthrosis of the shoulder	54
Arthrosis of the ankle joints	55
Baker's cysts	56
Prepatellar bursitis	57

Chondromalacia patellae, Chondropathia patellae	58
Distorsions in the hip region	59
Radial epicondylitis	60
Fractures in the hip region	61
Contusions in the hip region	62
Bekhterev's arthritis	63
Muscular tensions, , cervical and upper thoracic spine	64
Myogeloses	65
Osteoporosis	66
Syndrome of the apex of patella	67
Periarthritis humeroscapularis	68
Postoperative intestinal atonies	69
Radicular syndromes	70
Pain in the iliosacral joints	71
SUDECK's syndrome	72
Acute torticollis	73
Psychiatry	74
Anxiety disorder , Depressions	74
Withdrawal symptoms in alcoholism Withdrawal symptoms in drug addiction	75
Phobiae	76
Index	77

Origin of the HighTone Power Therapy

HighTone power therapy is a new quantum leap in the field of electro therapy.

Conventional electro therapy stimulates nerves and muscles. The primary aim of HighTone power therapy is a direct effect on the metabolism of the cells.

The differences can be seen in **effect**, **technology** and **practical use**.

The differences in effect relate to two main mechanisms of HighTone power therapy: Firstly, it introduces energy into the body to increase the cell energy potential. Secondly, it makes the cell structures oscillate to normalise the metabolism.

From the technical point of view the classical electro therapy is a modulation of amplitude: The current intensity is modulated, but the frequency remains constant. Electro therapy uses modulation frequencies between 0 and 200 Hertz in the low frequency range, and mostly 4.000 Hertz in the medium frequency range as carrier frequency.

In HighTone power therapy the amplitude and the frequency are modulated simultaneously. The higher the frequency, the more energy can be introduced correlating to the individual threshold curve of the patient's electro-sensitivity.

Thus, it is a **SIMULTaneous Frequency and Amplitude Modulation**. In the following it will always be called **SimulFAM®i**.

The letter 'i' stands for intensity, which increases simultaneously with rising frequency.

In the HighTone power therapy device **HiToP®** uses HighTone frequencies between 4.096 and 32.768 Hertz.


■ HighTone power therapy differs in two further practical points:

1. Classical electro therapy is applied for 5 to 10 minutes, but HighTone power therapy up to 60 minutes.
2. In classical electro therapy, as a rule, one channel and two electrodes are used. In HighTone power therapy up to 4 channels are in parallel operation, with 10 electrodes.

What are the two main effects of HighTone Power Therapy?

■ Introduction of energy into the body in order to activate the cells and to vitalise the body.

■ Creation of resonances in order to make the cells and the tissue structures oscillate, with the aim to improve metabolic processes and to relieve pain.



HiToP® HighTone Power Therapy

In 1935, Siegfried Koeppen already thought about the possibilities to use tone frequency therapy.

Dr. med. Hans-Ulrich May is recognized as the „father“ of the modern HighTone Therapy. Since 1988, he has been studying the variously effective applications very successfully.

Therapeutic efficiency

Already after the first treatment, 60 – 80 % of the patients do not have any troubles anymore or experience a great relief of pain (according to a patients' survey in 2001).

Practical efficiency

Patients with various kinds of pain syndromes (e.g. gonarthrosis, coxarthrosis, lumbar spine syn-drome, cervical spine syndrome, shoulder-arm-syndrome, etc.) experience significant relief of their pains which could not be achieved by conventional therapy methods by now.

This also means fast mobilization as well as the return to the required powerfulness for the working patient who has a job. Elderly people will feel much more vitality and their quality of life will increase.

Resonance phenomena


Apply an electric alternating field

Charged particles (molecules, cell structures, ...) of different sizes without an electric field (original position)




Resonance: Rising Frequency

32 kHz



Resonance: Sinking Frequency

8 kHz



What makes the cell and tissue structures oscillate?

The applied frequencies range from 4.096 to 32.768 Hertz over 3 octaves in 72 quarter-tone steps of 1 second each. The tones in the three octaves are the same that are used in standard music.

These high tone frequencies between 4.096 and 32.768 Hertz pass through the body in form of an electrical field that makes the charged particles oscillate.

The frequencies of the oscillations introduced create resonances in the molecules and cell structures.

Different frequencies activate structures of different size. For this reason it is important to offer a broad spectrum of frequencies.

The oscillations of the different particles in the tissue lead to many effects. One of them is a strongly increased distribution of pain and inflammation mediators as well as a positive effect on the transport of nutritive and waste substances.

Thus, the result is an improvement of metabolic processes, purging and pain relief.

How is the energy introduced into the cells?

The threshold curve of the electrosensitivity defines how much energy can be introduced into the body. The higher the frequency of the energy entered, the more energy the human body is able to tolerate.

The goal is to transfer as much energy as possible to the body. This can be achieved with high frequencies. The HighTone power therapy operates in the range between 4.096 and 32.768 Hertz. Up to 5.000 mWatt are introduced directly into the body.

HighTone power therapy increases the mitochondria in amount and size, the so-called energy „power stations“.

After the therapy the patient feels refreshed, just like after a long walk.

This energy input is made twice: once locally in the area to be treated and simultaneously on the whole body.

The whole-body vitalisation is important, for a disease is rarely limited to one area. The resistance of the whole body is improved.

HiToP[®] programs

The programs of the **HiToP** are pre-programmed sequences of treatments that are designed to achieve well defined effects for the customer.

In HiToP[®] 184The user can chose between the Program menu or the Indication menu by selecting the corresponding menu in the device settings.

The different programs can typically be used with one application of the electrodes, i.e. it is not necessary to change electrode application as the program changes from one treatment to the next.

The following is a brief description of the different programs and their typical use.

Vitalization

Vitalization consists of SimulFAM[®] *i* only and is typically applied as the whole-body treatment.

The program length is pre-set to 30 minutes but can be changed by the user.

Lipolysis

Lipolysis is the body's own process of freeing fatty cells from the tissue and making them accessible as energy carriers for the body to burn. This Lipolysis can be stimulated by HighTone Therapy. The scientifically proven frequency for maximum stimulation of lipolysis is 3 Hz.

In order to avoid adaptation of the body and based on practical experience our Lipolysis program uses a frequency sweep from 2.5 – 3.5 Hz. The length of the treatment is preset to 30 minutes but can be changed by the user.

As this program will only free but not burn the freed fatty cells it is recommended to follow it up with physical exercise.

The program length is pre-set to 30 minutes but can be changed by the user.

Cellulite

The cellulite treatment combines 10 minutes of lipolysis with 5 minutes of a frequency sweep from 10 – 20 Hz. This sweep will result in a vasoconstriction (10 Hz), pressing blood out of the area, and strong muscle stimulation (20 Hz), burning the freed fatty cells.

The program is of a fixed length of 30 minutes and the above described sequence is carried out twice.

Fatburning

This program focuses on the reduction of fatty cells. Therefor 80% of the treatment time is used for Lipolysis. The rest of the time a frequency sweep of 10 – 30 Hz is applied, giving some vasoconstriction but mainly muscle work in order to burn the freed fatty cells.

The program is of fixed length 30 minutes of 3 cycles with 8 minutes Lipolysis and 2 minutes 10 – 30 Hz.

Body Shaping

Goal of this program is to bring the body into a good shape by reduction of fatty cells on one side and muscle build up on the other side. Therefor we have evenly combined Lipolysis treatment with a efficient training of the muscles with SimulFAM[®] *X* frequencies between 20 and 30 Hz.

Important for a effective muscle training is a strong contraction. To achieve this the intensity is increased during the muscle stimulation period. To give the muscle time to relax, pauses are inserted corresponding to the Sport program.

The program is of fixed length 30 minutes of 2 cycles. Each cycle contains 5 minutes Lipolysis, 5 minutes muscle work and 5 minutes relaxation of the muscles with 150 – 200 Hz.

Electro-Massage

To achieve a relaxing and comforting loosening of a muscle or muscle group we can apply a frequency sweep of 0.5 to 3 Hz. This corresponds to the rhythm of a manual massage and is therefor called Electro-massage.

The program length is pre-set to 30 minutes but can be changed by the user.

Shaking-Massage

To achieve a relaxing and comforting loosening of a muscle or muscle group we can apply a frequency sweep of 3 - 10 Hz. This will result in a vibration or shaking of the muscle and is therefor called Shaking-massage.

The program length is preset to 30 minutes but can be changed by the user.

Sport

This is strong muscle training pure. The muscle is being brought to a strong tetanic contraction using SimulFAM® X frequencies between 20 and 30 Hz. The intensity should be clearly above the motoric threshold to ensure a strong contraction.

In order to let the muscle have the necessary relaxation pauses the training follows the cycle 5 s rise time, 5 s holding time (contraction) and 20 s pause.

A important innovation that gbo HighTone Therapy adds to the muscle training is that we can use a sub-threshold 32 kHz current (i.e. the user does not feel it) during the pauses to alleviate the metabolism and thereby greatly improve the muscles ability to be supplied with new nutrients and to get rid of waste products. This allows us to shorten the pause times and to achieve a much more efficient training effect.

The length of the program is pre-set to 15 minutes but can be changed by the user.

Reha

The Reha program is designed for patients in rehabilitation that need a slow re-introduction to joint movements. Therefor we stimulate the corresponding muscles with frequencies varying slowly between 0.2 and 0.5 Hz.

With sufficient intensity this will generate physiological contraction yielding a slow and energy saving limb movement.

The program length is pre-set to 30 minutes but can be changed by the user.

Mobilization

The program Mobilization addresses the patient after the rehabilitation or in an advanced status where the motions should slowly be increased in speed. The program therefor consists of 14 s slow motion with 0.5 Hz followed by an increase in speed with a frequency going from 1 to 8 Hz in 21 seconds.

This cycle of 14s + 21s = 35s is then repeated until the end of the treatment time.

The program length is pre-set to 30 minutes but can be changed by the user.

Muscle Pump

Muscle pump is a slow muscle movement where the repetitive muscle contractions serve as a “muscle pump” to squeeze liquids out of the treated area. This is generated using 0,3 Hz constantly.

The program length is pre-set to 30 minutes but can be changed by the user.

Trauma

To treat patients directly after injuries we have combined the edema treatment followed by 100Hz for pain treatment (counter irritation) and muscle relaxation.

The program length is pre-set to 30 minutes but can be changed by the user.

Lymph-Drainage

This program uses muscle pump (0.3 Hz), vasoconstriction (10 Hz) and a slight muscle contraction (10 – 30 Hz) of 10 minutes each in sequence to achieve an efficient lymphatic drainage.

The program length is fixed at 30 minutes.

Edema

The Edema reduction program consists simply of SimulFAM® X 10 Hz, being the optimum frequency for vasoconstriction in order to reduce edemas efficiently.

The program length is pre-set to 30 minutes but can be changed by the user.

Pain

To relieve the patient of pain the unit applies a frequency sweep of 90 – 110 Hz using the principle of counter irritation to suppress the pain sensation. It is always applied between the point of pain and the CNS.

The program length is pre-set to 30 minutes but can be changed by the user.

Nerveblock

The pain block uses a fixed frequency of 12 kHz. It is being applied directly to the pain point where it is fatiguing the pain receptors (nociceptors).

Universal

This program is a scan through all frequencies from 0.1 to 200 Hz.


The program length is pre-set to 30 minutes but can be changed by the user.

Diab. PNP

This program is used for pain relief at diabetic polyneuropathies. Caused by muscle activation, the metabolism is improved, thus improving the sensibility to insulin.

The program length is pre-set to 60 minutes but can be changed by the user.

Neurodermatitis

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel II				
Channel III				
Channel IV				


Position: lying

Therapies:

Number 10
Frequency 2-3 times the week

Remarks:

Psoriasis

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel II				
Channel III				
Channel IV				

Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Radiation injuries

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, big		
	Place two or more electrodes around the lesion of the skin..		
Channel III			
Channel IV			



Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Adnexitis, chronic

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Two dorsally in the sacral region, two ventrally in the lower abdominal region.		
Channel III			
Channel IV			



Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Cystitis, Cystopyelitis, Cystopyelonephritis, chronische

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Eine bis zwei abdominal gegen eins bis drei dorsal (Kreuzbein, LWS, untere BWS)		
Channel III			
Channel IV			



Position: lying

Therapies:


Number 10

Frequency 2-3 times a week

Remarks:

Dysmenorrhea

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Two dorsally in the sacral region, two ventrally in the lower abdominal region.		
Channel III			
Channel IV			



Position: lying

Therapies:


Number 10

Frequency 2-3 times a week

Remarks:

Pain caused by IUPs

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	One dorsally against one abdominally		
Channel III			
Channel IV			



Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Stress-urine incontinence

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	One ventrally abdominally two on the backside		
Channel III			
Channel IV			







Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Laryngitis , chronic

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment			
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel II	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, self adhesive electrodes			
	A larger one onto the cervix, a smaller one onto the larynx.			
Channel III				
Channel IV				



Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Chronic otitis media

Channel I	SimuFAM® i	60 min	 
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, medium size with spong pockets		
	One electrode on each ear. Use with well moistened spong pockets.		
Channel II			
Channel III			
Channel IV			


Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Frontal sinusitis , chronic

Channel I	SimuFAM® i	60 min	
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes or self adhesive electrodes		
	One electrode frontally above the nasal root one electrode suboccipitally.		
Channel II			
Channel III			
Channel IV			



Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Maxillary sinusitis , chronic

Channel I	SimuFAM® i	60 min	 
	<i>Settings at 4 and 16 kHz at threshold</i>		
	Plate electrodes or self adhesive electrodes		
	One electrode each to the right and left infraorbitally.		
Channel II			
Channel III			
Channel IV			


Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Tinnitus

Channel I	SimuFAM® i	60 min	
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, medium		
	Binaurally (one each on both auricles). Use with well moistened spong pockets.		
Channel II			
Channel III			
Channel IV			


Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Adipositas, overweight

Channel I	SimuFAM® X 3Hz/Lipolysis		60 min	
	Settings motorically above threshold			
	Plate electrodes, large			
	One pair each to the right and left of the waist			
Channel II				
Channel III				
Channel IV				

Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Arterial hypertonia

Arterial occlusive diseases

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel II				
Channel III				
Channel IV				

Position: lying


Therapies:

Number 10

Frequency 2-3 times a week

Remarks:

Uratic arthritis

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Two electrodes under the soles of foot, two electrodes on the palms.			
Channel II				
Channel III				
Channel IV				

Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Spasmodic asthma

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Two each cranially between scapula and cervical spine, two ventrally in the thoracoabdominal transitional region.		
Channel III			
Channel IV			



Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Decubitus in the back region

Pressure gangrene

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment			
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel II	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Around the pressure gangrene or the ulcer			
Channel III				
Channel IV				

Position: lying



Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Decubitus in the region of the heel bone

Pressure gangrene

Channel I	SimuFAM® i	60 min	
	<i>Settings at 4 and 16 kHz at threshold</i>		
	Plate electrodes, large		
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II			
Channel III			
Channel IV			


Position: lying

Therapies:

Number 10
 Frequency 2-3 times a week

Remarks:

Diabetic angiopathy

Channel I	SimuFAM® i	60 min	
	<i>Settings at 4 and 16 kHz at threshold</i>		
	Plate electrodes, large		
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II			
Channel III			
Channel IV			

Position: lying


Therapies:

Number 10
 Frequency 2-3 times a week

Remarks:

Thromboangiitis obliterans

Channel I	SimuFAM® <i>i</i>		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® X 100 Hz/pain		60 min
	Seettings slightly above threshold		
	Plate electrodes, large		
	One up to two lumbally against two plantarly		
Channel III			
Channel IV			





Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Hematoma in the sural region

Channel I	SimuFAM® i	60 min	
	<i>Settings at 4 and 16 kHz at threshold</i>		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® X 10 Hz/Edema	60 min	
	<i>Settings slightly above threshold</i>		
	Plate electrodes, large		
	One electrode at calf and one at cnemis		
Channel III			
Channel IV			



Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Cardiac discomforts, functional

Channel I	SimuFAM® i	60 min	
	<i>Settings at 4 and 16 kHz at threshold</i>		
	Plate electrodes, large		
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II			
Channel III			
Channel IV			



Position: lying

Therapies:

Number 10
 Frequency 2-3 times a week

Remarks:

Migraine (treatment at intervals)

Channel I	SimuFAM® i	60 min	
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes or self adhesive electrodes		
	Place one electrode or special four-button electrode in the upper neck suboccipitally, the counter-electrode caudally to the first one in the lower cervical spine region		
Channel II			
Channel III			
Channel IV			



Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Migraine and BING-HORTON syndrome

Channel I	SimuFAM® X 100 Hz/pain	max. 30 min	
	<i>Setting on threshold</i>		
	Plate electrode or four-button electrode		
	One electrode or special four-button electrode in the region of the carotid sinus		
	do not fix, but have it held! One counter-electrode nuchally..		
Channel II			
Channel III			
Channel IV			

Position: lying

Therapies:


Number 10
 Frequency 2-3 times a week

Remarks:

Obstipation atonic

Obstipation, spastic

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Two up to four abdominally against one up to two dorsally (sacrum, lumbar spine, lower thoracic spine)		
Channel III			
Channel IV			



Position: lying



Therapies:

Number 10

Frequency 2-3 times a week

Remarks:

Posttraumatic edema in the area of the hand

Channel I	SimuFAM® X 10 Hz/Edema	60 min	
	<i>Settings slightly above threshold</i>		
	Plate electrodes, large		
	Proximal treatment: One electrode medially (axially) in the area of the vessel-nerve cord one counter-electrode laterally (radially).		
Channel II			
Channel III			
Channel IV			


Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Posttraumatic edema in the area of the foot

Channel I	SimuFAM® X 10 Hz/Edema	60 min	
	Settings slightly above threshold		
	Plate electrodes, large		
	Proximal treatment: One electrode medially (axially) in the area of the vessel-nerve cord one counter-electrode laterally (radially).		
Channel II			
Channel III			
Channel IV			



Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Raynaud's syndrome

Channel I	SimuFAM® X 100 Hz/pain	60 min	
	<i>Settings slightly above threshold</i>		
	Plate electrodes, large		
	One large electrode in the neck against two palmar electrodes placed on both sides		
Channel II			
Channel III			
Channel IV			

Position: lying





Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Thrombosis prophylaxis

With HiToP 184 you can do a whole body treatment for vitalisation at the same time.(Channel III)

Channel I	SimuFAM® X 0,3 Hz/Musclepump		60 min	
	Settings motorically above threshold			
	Plate electrodes or self adhesive electrodes			
	Place two electrodes each on the crural regions proximally and distally, laterally to the anterior margin of the tibia onto the foot and toe muscle extensor			
Channel II	SimuFAM® X 0,3 Hz/Musclepump		60 min	
	Settings motorically above threshold			
	Plate electrodes or self adhesive electrodes			
	Place two electrodes each on the crural regions proximally and distally, laterally to the anterior margin of the tibia onto the foot and toe muscle extensor			
Channel III	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel IV				

Position: lying


Therapies:

Number 10

Frequency 2-3 times a week

Remarks:

Stasis ulcer

Channel I	SimuFAM® i	60 min	
	<i>Settings at 4 and 16 kHz at threshold</i>		
	Plate electrodes, large		
	Two electrodes under the soles of foot.		
Channel II			
Channel III			
Channel IV			

Position: lying





Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Venous insufficiency , chronic

With HiToP 184 you can do a whole body treatment for vitalisation at the same time.(Channel III)

Channel I	SimuFAM® X 0,3 Hz/Musclepump 60 min		
	Settings motorically above threshold		
	Plate electrodes or self adhesive electrodes		
	Place two electrodes each on the crural regions proximally and distally, laterally to the anterior margin of the tibia onto the foot and toe muscle extensor		
Channel II	SimuFAM® X 0,3 Hz/Musclepump 60 min		
	Settings motorically above threshold		
	Plate electrodes or self adhesive electrodes		
	Place two electrodes each on the crural regions proximally and distally, laterally to the anterior margin of the tibia onto the foot and toe muscle extensor		
Channel III	SimuFAM® i	60 min	
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel IV			


Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Amyotrophic lateral sklerosis

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel II				
Channel III				
Channel IV				


Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Atypical faciocephalgia

Channel I	SimuFAM® X 100 Hz/Schmerz	60 min	
	<i>Settings slightly above threshold</i>		
	Plate electrodes, large		
	Two electrodes under the soles of foot for counter irritation		
Channel II			
Channel III			
Channel IV			


Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Diabetic polyneuropathy

Channel I	SimuFAM® X 20 Hz/Diab. PNP	60 min	
	<i>Settings motorically above threshold</i>		
	Plate electrodes 120 x 80		
	Each one electrode to both thighs to the knee, each one electrode to both thighs to the grion		
Channel II			
Channel III			
Channel IV			

Position: lying

Therapies:


Number 10
Frequency 2-3 times a week

Remarks:

Nocturnal enuresis

bedwetting

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	One ventrally abdominally, two in the proximal medio-dorsal upper thigh regions below the gluteal folds..		
Channel III			
Channel IV			




Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Facial nerve paresis

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes or self adhesive electrodes			
	One flexible electrode on the affected half of the face, a larger electrode in the nuchal region (dorsal cervical spine region).			
Channel II				
Channel III				
Channel IV				


Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Polyneuropathies

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel II				
Channel III				
Channel IV				


Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Shoulder pain in hemiplegia

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	One ventrally, one dorsally to the shoulder joint			
Channel II	SimuFAM® X 100 Hz/pain		60 min	
	Settings slightly above threshold threshold			
	Plate electrodes			
	Two electrodes on the palm.			
Channel III				
Channel IV				


Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Tension headache, Vasomotoric headache

Channel I	SimuFAM® X 100 Hz/pain		60 min	
	Settings slightly above threshold threshold			
	Plate electrodes or four button electrode			
	Place one electrode or special four-button electrode to the upper neck suboccipitally, the counter-electrode caudally to the first one in the lower cervical spine region or the upper thoracic spine region			
Channel II				
Channel III				
Channel IV				

Position: lying



Therapies:

Number 10

Frequency 2-3 times a week

Remarks:

Conjunctivitis, allergic

Channel I	SimuFAM® X 10 Hz/Ödem	60 min	
	<i>Settings at threshold</i>		
	Plate electrode, small or self adhesive electrodes		
	One electrode each onto the eyes.		
Channel II			
Channel III			
Channel IV			



Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Arthrosis of the carpus

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment			
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel II	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes or self adhesive electrodes			
	One electrode volarly, one electrode dorsally			
Channel III				
Channel IV				

Position: lying




Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Arthrosis of the hip joints, Coxarthrosis

With HiToP 184 you can do a whole body treatment for vitalisation at the same time.(Channel III)

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Two electrodes, one each to the right and left of the hip joint region..			
Channel II	SimuFAM® X 100 Hz/pain		60 min	
	Settings slightly above threshold			
	Plate electrodes, large			
	Two or four electrodes in unilateral or bilateral treatment: One each proximally dorsally lumbally One each distally ventrally on the upper thigh..			
Channel III	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel IV				

Position: lying





Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Arthrosis of the knee joints

With HiToP 184 you can do a whole body treatment for vitalisation at the same time.(Channel III)

Channel I	SimuFAM® i	60 min	
	Settings slightly above threshold		
	Plate electrodes, large		
	One electrode on both sides of the knee joint.		
Channel II	SimuFAM® X 100 Hz/pain	60 min	
	Einstellung slightly überschwellig		
	Plate electrodes, large		
	One electrode on the thigh, one electrode on the calf.		
Channel III	SimuFAM® i	60 min	
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel IV			

Position: lying





Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Arthrosis of the shoulder

With HiToP 184 you can do a whole body treatment for vitalisation at the same time.(Channel III)

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes or self adhesive electrodes			
	One electrode ventrally, one electrode dorsally to the shoulder joint.			
Channel II	SimuFAM® X 100 Hz/pain		60 min	
	Settings slightly above threshold			
	Plattenelektroden möglichst groß			
	One electrodes on each upper arm.			
Channel III	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel IV				

Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Arthrosis of the ankle joints

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, medium		
	One tibially, one fibularly.		
	Use with well moistened spong pockets.		
Channel III			
Channel IV			



Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Baker's cysts

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	One electrode on the knee cap, one electrode in the hollow of the knee.		
Channel III			
Channel IV			



Position: lying



Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Prepatellar bursitis

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	One electrode on the knee cap, one electrode in the hollow of the knee.		
Channel III			
Channel IV			







Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Chondromalacia patellae, Chondropathia patellae

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel II	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	One electrode on the knee cap, one electrode in the hollow of the knee.			
Channel III				
Channel IV				

Position: lying

Therapies:


Number 10

Frequency 2-3 times a week

Remarks:

Distorsions in the hip region

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® X 10 Hz/Trauma		60 min
	Settings slightly above threshold		
	Plate electrodes, large		
	Two electrodes, one each to the right and left of the hip joint region.		
Channel III			
Channel IV			



Position: lying





Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Radial epicondylitis

With HiToP 184 you can do a whole body treatment for vitalisation at the same time.(Channel III)

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes			
	One (smaller) electrode radially onto the epicondylus, one (larger) or two counter-electrodes ulnarly to the opposite			
Channel II	SimuFAM® X 100 Hz/pain		60 min	
	Settings slightly above threshold			
	Plate electrodes, large			
	One electrode tu the upper arm and one electrode th the forearm			
Channel III	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel IV				

Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Fractures in the hip region

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® X 10 Hz/Ede ma		60 min
	Settings slightly above threshold		
	Plate electrodes, large		
	Two electrodes, one each to the right and left of the hip joint region..		
Channel III			
Channel IV			



Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Contusions in the hip region

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® X 10 Hz/Trauma		60 min
	Settings slightly above threshold		
	Plate electrodes, large		
	Two electrodes, one each to the right and left of the hip joint region...		
Channel III			
Channel IV			





Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Bekhterev's arthritis

Channel I	SimuFAM® i	60 min	
	<i>Settings at 4 and 16 kHz at threshold</i>		
	Plate electrodes, large		
	Four electrodes with changing polarity uniformly distributed to the spine area.		
Channel II			
Channel III			
Channel IV			


Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Muscular tensions, , cervical and upper thoracic spine

Channel I	SimuFAM® X 1,5 – 10 Hz /Shaking	60 min	
	massage		
	Settings motorically above threshold		
	Plate electrodes or self adhesive electrodes		
Channel II	One plate electrode or self adhesive electrode each onto the upper edge of the trapezius muscle on both sides		
Channel III			
Channel IV			

Position: lying



Therapies:

Number 10

Frequency 2-3 times a week

Remarks:

Myogeloses

Channel I	SimuFAM® i	60 min	
	<i>Settings at 4 and 16 kHz at threshold</i>		
	Plate electrodes, large		
	One plate electrode or self adhesive electrode each onto the upper edge of the trapezius muscle on both sides		
Channel II			
Channel III			
Channel IV			

Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Osteoporosis

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Two ventrally abdominally against two three or four dorsally along the spine from cervical spine upto the sacrum..		
Channel III			
Channel IV			



Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Syndrome of the apex of patella

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	One on the knee cap, one on the pople.		
Channel III			
Channel IV			



Position: lying





Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Periarthritis humeroscapularis

With HiToP 184 you can do a whole body treatment for vitalisation at the same time.(Channel III)

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes or self adhesive electrodes			
	One electrode ventrally, one electrode dorsally to the shoulder joint.			
Channel II	SimuFAM® X 100 Hz/pain		60 min	
	Settings slightly above threshold			
	Plate electrodes, large			
	One electrode to both upper arms..			
Channel III	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel IV				

Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Postoperative intestinal atonies

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	One to two abdominally against one to two dorsally (sacrum, lumbar spine, lower thoracic spine)		
Channel III			
Channel IV			



Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Radicular syndromes

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Two to four electrodes ventrally abdominally against two electrodes dorsally in the area of the lower thoracic and the lumbar spine		
Channel III			
Channel IV			



Position: lying


Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Pain in the iliosacral joints

Channel I	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment		
	Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		60 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	One electrode each on both ileosacral joints		
Channel III			
Channel IV			




Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

SUDECK's syndrome

Channel I	SimuFAM® i		60 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel II				
Channel III				
Channel IV				


Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Acute torticollis

Channel I	SimuFAM® X 100 Hz/Schmerz	60 min			
	<i>Settings slightly above threshold</i>				
	Plate electrodes or self adhesive electrodes				
	Place one transverse suboccipitally one transverse below that (caudal direction) in the neck				
Channel II					
Channel III					
Channel IV					


Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Anxiety disorder , Depressions

Channel I	SimuFAM® i		30 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel II	SimuFAM® i		30 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes			
	Binaurally (one each on both auricles). Use electrodes with good moistened spong pockets.			
Channel III				
Channel IV				

Position: lying

Therapies:


Number 10
Frequency 2-3 times a week

Remarks:

Withdrawal symptoms in alcoholism

Withdrawal symptoms in drug addiction

Channel I	SimuFAM® i		30 min
	Settings at 4 and 16 kHz at threshold		
	Plate electrodes, large		
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck		
Channel II	SimuFAM® i		30 min
	Einstellung bei 4 und 16 kHz schwellig		
	Plate electrodes medium with spong pockets		
	Binaurally (one each on both auricles). Use electrodes with good moistened spong pockets.		
Channel III			
Channel IV			



Position: lying


Therapies:

Number 10

Frequency 2-3 times a week

Remarks:

Phobiae

Channel I	SimuFAM® i		30 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes, large			
	Whole-body treatment Two electrodes under the soles of foot, two electrodes on the forearms, one electrodes in the neck			
Channel II	SimuFAM® i		30 min	
	Settings at 4 and 16 kHz at threshold			
	Plate electrodes medium with spong pockets			
	Binaurally (one each on both auricles). Use electrodes with good moistened spong pockets.			
Channel III				
Channel IV				

Position: lying

Therapies:

Number 10
Frequency 2-3 times a week

Remarks:

Index

- Acute torticollis 72
- Adipositas 23
- Adnexitis 13
- Amyotrophic lateral sklerosis 42
- Anxiety disorder 73
- Arterial hypertonia 24
- Arterial occlusive diseases 24
- Arthrosis of the carpus 50
- Arthrosis of the ankle joints 54
- Arthrosis of the hip joints 51
- Arthrosis of the knee joints 52
- Arthrosis of the shoulder 53
- Atypical faciocephalalgia 43
- Baker's cysts 55
- Bekhterev's arthritis 62
- BING-HORTON syndrome 34
- Cardiac discomforts, functional 32
- Chondromalacia patellae 57
- Chondropathia patellae 57
- Chronic otitis media 19
- Conjunctivitis, allergic 49
- Contusions in the hip region 61
- Coxarthrosis 51
- Cystitis 14
- Cystopyelitis 14
- Cystopyelonephritis 14
- Decubitus 27
- Depressions 73
- Diabetic angiopathy 29
- Distorsions in the hip region 58
- Dysmenorrhea 15
- epicondylitis 59
- Facial nerve paresis 45
- Fractures in the hip region 60
- Frontal sinusitis 20
- Hematoma in the sural region 31
- Laryngitis 18
- Maxillary sinusitis 21
- Migraine 34
- Migraine (treatment at intervals) 33
- Muscular tensions 63
- Myogeloses 64
- Neurodermatitis 10
- Nocturnal enuresis 44
- Obstipation 35
- Osteoporosis 65
- Pain caused by IUPs 16
- Pain in the ileosacral joints 70
- Periarthritis humeroscapularis 67
- Phobias 75
- Polyneuropathies 46
- Postoperative intestinal atonies 68
- Posttraumatic edema 36, 37
- Prepatellarbursitis 56
- Psoriasis 11
- Radiation injuries 12
- Radicular syndromes 69
- Raynaud's syndrome 38
- Shoulder pain in hemiplegia 47
- Spasmodic asthma 26
- Stasis ulcer 40
- Stress-urine incontinence 17
- SUDECK's syndrome 71
- Syndrome of the apex of patella 66
- Tension headache 48
- Thromboangiitis obliterans 30
- Thromboseprophylaxe 39
- Tinnitus 22
- Uratric arthritis 25
- Vasomotoric headache 48
- Venous insufficiency 41
- Withdrawal symptoms in alcoholism 74
- Withdrawal symptoms in drug addiction 74

Notes

[illegible]

Notes

[illegible]

011-7-0116

gbo Medizintechnik AG
Kleiststrasse 6
D-64668 Rimbach
Germany

Tel.: +49(0)6253/808-0
Fax: +49(0)6253/808-245
e-mail: Info@gbo-med.de
internet: www.gbo-med.de