

TECHNICAL INFORMATION NOTICE

Following products have been tested for event-related sterility maintenance study, according to the use they are destined to:

- TYVEK®/FILM BAG for hydrogen peroxide sterilization
- TYVEK®/FILM ROLL for hydrogen peroxide sterilization

Test results after 180 days are correct and samples satisfy sterility test.

This information is based on the results given by BIOTECNAL laboratory which has performed the test according to procedures of UNE-EN-ISO 11737-2.

Greetings

Quality Control



Sabadell, 10th May 2010