

# medical packaging



# medical packaging

Excellent quality packaging to fulfil the most exigent sterilization processes. Designed to be used in Hospitals, Clinics, Sterilization Centres, and Vet Centres as well as in medical environments.

Manufactured and labelled according to ISO11607-1 (EN-868), including barcode in product labels.

At the same time, manufacturing process and service are according to procedures of ISO 9001-2008.

Design and Quality mix in this product, providing an excellent peelability and easy opening. Film transparency shows product without needing of opening and confers an attractive appearance to the product inside.

Printing in several languages, makes it easy to deal and sell in several markets and for several users.

All our products are in stock and ready to be expedited worldwide in 72 hours.

Special sizes can be manufactured, upon request.

○ service

○ availability

○ design

○ peelability

○ quality

○ versatility

# packaging for steam/eo/formol

Pouches and rolls, manufactured with 60 g/m<sup>2</sup> medical paper and a plastic film (blue coloured). Indicators for Steam/EO/Formol are also printed in packaging.

## paper/film flat bag

Code	Sizes								U/Box	Box cm	Box/KG
MBO1000	75	x	210	mm	3,0	"	8,3	"	2.500	31x29x11	4,8
MBO1001	75	x	280	mm	3,0	"	11,0	"	2.400	31x29x11	6,5
MBO1003	100	x	200	mm	3,9	"	7,9	"	2.400	41x29x11	5,9
MBO1004	100	x	240	mm	3,9	"	9,4	"	2.400	41x29x11	7
MBO1005	100	x	280	mm	3,9	"	11,0	"	2.400	41x29x11	8,5
MBO1006	100	x	420	mm	3,9	"	16,5	"	2.400	41x43x11	12,5
MBO1020	100	x	550	mm	3,9	"	21,7	"	2.000	41x43x11	13,2
MBO1022	100	x	580	mm	3,9	"	22,8	"	2.000	41x59x11	14,2
MBO1007	150	x	210	mm	5,9	"	8,3	"	2.000	32x43x11	8
MBO1008	150	x	280	mm	5,9	"	11,0	"	2.000	31x57x11	10,5
MBO1009	150	x	320	mm	5,9	"	12,6	"	2.000	37x65x11	11,5
MBO1010	150	x	360	mm	5,9	"	14,2	"	2.000	37x65x11	13
MBO1021	150	x	400	mm	5,9	"	15,7	"	1.200	32x43x11	8,7
MBO1011	200	x	320	mm	7,9	"	12,6	"	1.800	37x65x11	14
MBO1012	200	x	360	mm	7,9	"	14,2	"	1.800	37x65x11	15,8
MBO1013	200	x	420	mm	7,9	"	16,5	"	1.200	41x43x11	12,5
MBO1023	200	X	500	mm	7,9	"	19,7	"	1.200	51x51x11	14,8
MBO1014	250	x	360	mm	9,8	"	14,2	"	1.200	51x43x11	13,5
MBO1015	250	x	420	mm	9,8	"	16,5	"	1.200	51x43x11	15,5
MBO1016	250	x	500	mm	9,8	"	19,7	"	1.200	51x51x11	18,5
MBO1017	300	x	420	mm	11,8	"	16,5	"	600	32x43x11	9,5
MBO1018	300	x	500	mm	11,8	"	19,7	"	600	31x57x11	11
MBO1024	300	x	580	mm	11,8	"	22,8	"	700	41x59x11	14,2
MBO1019	400	x	580	mm	15,7	"	22,8	"	600	41x59x11	17



paper/film  
flat bag



paper/film  
gusset bag



paper/film  
flat roll

## paper/film gusset bag

Code	Sizes								U/Box	Box cm	Box/KG				
MBF1202	200	x	55	x	400	mm	7,9	x	2,2	"	15,7	"	600	43x26x23	7

paper/film selfsealing bag

Code	Sizes									U/Pack	U/Box	Box cm	Box/KG
MBA1080	75	x	200	mm	3,0	"	7,9	"		200	1.200	15x24x17	2,7
MBA1081	90	x	200	mm	3,5	"	7,9	"		200	1.200	28x27x12	4
MBA1082	90	x	230	mm	3,5	"	9,1	"		200	1.200	28x27x12	4,6
MBA1084	100	x	240	mm	3,9	"	9,4	"		200	1.200	31x29x11	4,5
MBA1085	130	x	250	mm	5,1	"	9,8	"		200	1.200	41x29x11	5,8
MBA1086	130	x	330	mm	5,1	"	13,0	"		200	1.200	41x37x11	7,6
MBA1087	140	x	280	mm	5,5	"	11,0	"		200	1.200	43x32x11	7
MBA1088	200	x	330	mm	7,9	"	13,0	"		200	800	41x37x11	7,8
MBA1089	250	x	400	mm	9,8	"	15,7	"		200	800	51x43x11	12
MBA1090	300	x	390	mm	11,8	"	15,4	"		***	400	32x43x11	6,9
MBA1091	300	x	450	mm	11,8	"	17,7	"		***	400	31x49x11	8



paper/film selfsealing bag



selfsealing tape



paper/film gusset roll

paper/film flat roll

Código	Sizes									U/Box	Box cm	Box/KG
MRP1040	50	mm	x	200	mts	2,0	"	656	Feet	6	31x25x23	6
MRP1041	75	mm	x	200	mts	3,0	"	656	Feet	4	31x25x23	8
MRP1042	100	mm	x	200	mts	3,9	"	656	Feet	4	43x26x23	11
MRP1043	125	mm	x	200	mts	4,9	"	656	Feet	4	43,5x32x23	14
MRP1044	150	mm	x	200	mts	5,9	"	656	Feet	4	43,5x32x23	17
MRP1045	200	mm	x	200	mts	7,9	"	656	Feet	2	43x26x23	11
MRP1046	250	mm	x	200	mts	9,8	"	656	Feet	2	52x28x23	14
MRP1047	300	mm	x	200	mts	11,8	"	656	Feet	2	43,5x32x23	17
MRP1048	400	mm	x	200	mts	15,7	"	656	Feet	1	43x26x23	11

paper/film gusset roll

Code	Sizes									U/Box	Box cm	Box/KG
MRF1060	100 x 40	mm	x	100	mts	3,9 x 1,6	"	328	Feet	4	43x26x23	7
MRF1061	150 x 40	mm	x	100	mts	5,9 x 1,6	"	328	Feet	4	43,5x32x23	11
MRF1062	200 x 60	mm	x	100	mts	7,8 x 2,3	"	328	Feet	2	43x26x23	7
MRF1063	250 x 60	mm	x	100	mts	9,8 x 2,3	"	328	Feet	2	52x28x23	8
MRF1064	300 x 80	mm	x	100	mts	11,8 x 3,1	"	328	Feet	1	31x25x23	6
MRF1065	400 x 80	mm	x	100	mts	15,7 x 3,1	"	328	Feet	1	43x26x23	7

# packaging for low temperature plasma (h<sub>2</sub> o<sub>2</sub>)/eo

Pouches and rolls, manufactured with Tyvek® and plastic film (transparent). Indicator for low temperature is printed in packaging.

## tyvek®/film flat bag

Code	Sizes								U/Box	Box cm	Box/KG
MBT1160	75	x	210	mm	3,0	"	8,3	"	1.000	16x22x15	2,4
MBT1161	100	x	280	mm	3,9	"	11,0	"	1.000	21x29x15	4,1
MBT1162	150	x	320	mm	5,9	"	12,6	"	1.000	31x33x15	7
MBT1163	200	x	420	mm	7,9	"	16,5	"	350	41x43x11	4
MBT1164	250	x	500	mm	9,8	"	19,7	"	350	51x51x11	6
MBT1165	320	x	560	mm	12,6	"	22,0	"	400	41x59x11	11,4
MBT1166	400	x	580	mm	15,7	"	22,8	"	350	41x59x11	11,5

## tyvek®/film flat roll

Code	Sizes								U/Box	Box cm	Box/KG	
MRT1170	50	mm	x	70	mts	2,0	"	230	Feet	6	31x25x23	3
MRT1171	75	mm	x	70	mts	3,0	"	230	Feet	4	31x25x23	3,5
MRT1172	100	mm	x	70	mts	3,9	"	230	Feet	4	43x26x23	4
MRT1174	150	mm	x	70	mts	5,9	"	230	Feet	4	43,5x32x23	6
MRT1175	200	mm	x	70	mts	7,9	"	230	Feet	2	41x21x22	4,8
MRT1176	250	mm	x	70	mts	9,8	"	230	Feet	2	52x28x23	5
MRT1177	300	mm	x	70	mts	11,8	"	230	Feet	2	43,5x32x23	6
MRT1178	400	mm	x	70	mts	15,7	"	230	Feet	1	43x26x23	4
MRT1179	500	mm	x	70	mts	19,7	"	230	Feet	1	53x26x23	5




tyvek®/film bag



opening instructions



tyvek®/film roll



plasma indicator tyvek®

## industrial packaging

Although our production line is focused in supplying stock products, we also print medical paper and films for manufacturers of medical devices for their packaging machines. Our experience and knowledge in manufacturing stock items, is an added value for our industrial customers, providing exactly what customer needs.


### products

- Medical paper.
- Pa/pp – pet/pp – pet/pe – pa/pe films.
- Tyvek®.
- Private label.

(Tyvek® is a registered trade mark of DUPONT)



BOLSAPLAST, S.L. C. Horta, 22. 08203 Sabadell (Barcelona) ESPAÑA  
Tel. +34 937117361 / Fax +34 937122355. Email: [medical@bolsaplast.com](mailto:medical@bolsaplast.com) / Website: [www.bolsaplast.com](http://www.bolsaplast.com)



steam/eo/formol  
packaging

low  
temperature  
plasma  
( $\text{h}_2 \text{o}_2$ )/eo  
packaging

industrial  
packaging

