

ELECTROTHERAPY OF THE NEXT GENERATION



# One Touch to Start

The new HiToP<sup>®</sup> 1touch

**gbo** 

# The new HiToP® 1touch Electrotherapy in a Class of its Own

## One Touch to Start - now also for mobile use

The HiToP® 1touch - the latest member of the innovative and unique HiToP® family of gbo Medizintechnik AG - provides local and holistic pain therapy as well as therapies for other indications. The user interface provides a maximum level of usability. With a single touch one can select a desired therapy out of the personal list of favorites. The therapy is started by simply turning the intensity regulator. The modern high-tech design made of brushed aluminium and the 10" TFT LCD full color touch screen monitor are an essential tool for the therapist. As with all other HiToP® devices the user can choose the individual language.

## Full Range of Functions and Ease of Use

The HiToP® 1touch is designed for the use in the doctor's practice as well as mobile at the patient's house. Via the touch screen the therapist can select either therapies or indications in a most easy fashion. The graphic display supports the work with the device: pictures of the application demonstrate the placement of the electrodes. The HiToP® 1touch can be used in two different positions: propped up and laid flat (see pictures). It can be easily and safely taken along with the matching transport bag. Thus the HiToP® 1touch is the ideal expansion to the HiToP® 4 and 2touch for the usage in the sports club or private environment.



## Technical Data

Treatment time	max. 90 minutes
Display	10.1" TFT
Power supply	100 – 240 Volt, 50 – 60 Hz
Power consumption	max. 45 VA
Class of protection	I
Class of equipment	IIa according to MDD
Dimensions	245 x 270 x 60 mm (h x w x d)
Weight	2.44 kg without accessories
Channels	1 channel
Warranty	24 months

Additional Information about the new HiToP® 1touch can be found on our homepage [www.gbo-med.de](http://www.gbo-med.de) as well as [www.hightonetherapy.com](http://www.hightonetherapy.com). Personal contacting is also welcome.

### gbo Medizintechnik AG

Kleiststrasse 6 · D-64668 Rimbach  
Telefon: +49 (0) 6253 / 808-0 · Fax: +49 (0) 6253 / 808-245  
E-Mail: [info@gbo-ag.de](mailto:info@gbo-ag.de) · [www.gbo-med.de](http://www.gbo-med.de)

