


ELECTROTHERAPY OF THE NEXT GENERATION



One Touch to Start

The new HiToP® Devices



The new HiToP® Devices

Electrotherapy in a Class of its Own

One Touch to Start

With its innovative and unique devices HiToP® 4touch and HiToP® 2touch gbo Medizintechnik AG presents to the user local and holistic therapies for different indications. Yet the user interface gives a maximum level of usability: With one single touch you can select a therapy from your personal list of favorites. Simply turning the intensity regulator allows a real quick-start. The modern high-tech design with brushed aluminum surface and the 15" TFT LCD full color touch screen monitor offer to the health care professional an indispensable feature.

Complete Range of Functions and Ease of Use

The therapeutic device comes with four or two independent channels which are being displayed individually. Via the touch screen the doctor can select the patient to be treated, therapies or indications in a most easy fashion. The graphic display supports the work with the device: pictures with applications demonstrate the placement of the electrodes. The device cart rounds out the usability concept, for each accessory has its defined position. A special highlight is the treatment bed which is especially designed for the treatment with HiToP® therapy and offers many details to render the treatment even more comfortable. Patients as well as health care professionals benefit from functions like quicker fixing of patient cables and electrodes, heated surface of the bed and the possibility to warm up the electrodes.



Technical Data

Treatment time	max. 90 minutes
Display	15" tft lcd
Power supply	100 – 240 Volt, 50 – 60 Hz
Power consumption	max. 200 Va
Class of protection	I
Class of equipment	Ila according to MDD
Dimensions	285 x 360 x 373 (H x W x D)
Weight	max. 13 kg
Independent channels	4 respectively 2
Warranty	24 months

You will find more detailed information on our latest HiToP® devices under www.gbo-med.de and under www.hightonetherapy.com. Personal contacting is also welcome.

gbo Medizintechnik AG

Kleiststrasse 6 · D-64668 Rimbach

Phone: +49 (0) 6253 / 808-0 · Fax: +49 (0) 6253 / 808-245

E-mail: info@gbo-ag.de · www.gbo-med.de

