



Prickle, burning, numbness?

HighTone Power Therapy for diabetics

HiToP®191

HIGHTONE POWER THERAPY

Finally free of pain

For patients suffering from polyneuropathy we offer our special unit HiToP®191 for treatment at home. Several clinical studies prove effectiveness and therapeutical efficiency of HighTone Power Therapy for diabetical and uremic polyneuropathy. The effects of this innovative therapy were evaluated in comparison to the conventional electrotherapy.

Under the direction of Prof. Dr. Stephan Martin the German Center of Diabetics (DDZ) Duesseldorf scientifically proved the efficiency of HighTone Power Therapy on diabetical polyneuropathy. In comparison to the conventional method of the Transcutaneous Electrical Nerve Stimulation (TENS), HighTone Power Therapy is much more effective. For 80 percent of the polyneuropatic patients treated with HighTone Power Therapy (painful and painless polyneuropathy) considerable relief was documented, only 33 percent of the patients treated with TENS experienced similar relief.

A clinical study under the responsibility of Prof. Dr. Peter P. Nawroth at the University Hospital in Heidelberg (Germany) investigated the effect of HighTone Power Therapy on patients with painful diabetical polyneuropathy: Also in this case considerable relief was documented for 65 percent of the patients.

An observational study by the West German Center of Diabetics and Health at Duesseldorf (WDGZ) extended the data base of the studies: 414 patients suffering from diabetes received a therapy unit for a 6 week treatment at home. For 88.4 percent a significant reduction of neuropatic pain was documented. Also sleep disorders were reduced.



Technical Data

display	LCD
output frequency	4096 – 32768 Hz
modulations	20 Hz
treatment time	1 – 90 minutes
dimensions	23 x 13,5 x 23 cm (W x H x D)
Weight	1.4 kg (without accessories)
Warranty period	24 months

you will find more detailed information under www.gbo-med.de .
of course you are always welcome to contact us personally.

gbo Medizintechnik AG

Kleiststrasse 6 · D-64668 Rimbach
Phone: +49 (0) 6253 / 808-0 · Fax: +49 (0) 6253 / 808-245
E-Mail: info@gbo-ag.de · www.gbo-med.de

